

Corporate Information

Your Partner...

Pace Medical... over 30 years of pacing experience.

Discover the Difference™

*...Solutions today for
Your pacing needs...*

*Technology
Reliability
Clinical Experience*



...designs and manufactures Temporary Cardiac Pacemakers, Pacing Analyzers and accessory products used in pacing the heart, diagnosing the performance of implantable pacing leads for permanent pacemakers and implantable defibrillators, while contributing to better patient care.

Corporate Profile...

Pace Medical is positioned to be the leading medical device company exclusively dedicated to designing and marketing a complete line of high-quality temporary cardiac pacing products to meet the expanding needs of the medical profession worldwide. Our products include both single and dual-chamber temporary cardiac pacemakers, a dual-chamber pacing analyzer, autoclavable and single-use temporary pacemaker patient extension and surgical cables and related accessories.

Products...

The AccuPace™, dual-chamber pacing analyzer, combines the benefits of dual-chamber, multi-mode, multi-parameter pacing with testing features that evaluate, display and store important characteristics of the patient's lead system. The temporary pacing feature allows the physician to perform pre-implant stimulation studies of the basic pacing parameters and functions on the patient prior to the implantation of a permanent pulse generator. Considerable attention has been paid to the "user-friendly" aspects of operation.

Our single-chamber temporary cardiac pacemakers include the Bedside™ model which attaches to the patient's bedrail and the Mini-Pacer™ which enables the patient to remain ambulatory during treatment. Both products are used worldwide and have significant clinical use.



AccuPace™ Pacing Analyzer



MICRO-PACE™ Dual-chamber Temporary Pacemaker



Bedside™ Single-chamber Temporary Pacemaker

Products...

The Company's MICRO-PACE™, dual-chamber temporary cardiac pacemaker, is software programmable and capable of multi-parameter, multi-mode operation. It has the ability to sense and pace at high rates in the DDD mode; thereby, restoring A-V synchrony and improving cardiac output.

The Pace Line™ series of autoclavable and single-use products include— surgical and patient extension cables and adapters— designed to enhance our ability to offer a complete line of highly-competitive temporary pacing products for our customers.

...Pace Medical, your answer to your pacing needs...



Universal™ Connector



Redel® Connector

P A C E M E D I C A L

Pace Medical, Inc.
391 Totten Pond Road
Waltham, MA 02451 USA
Tel: 781-890-5656
Fax: 781-890-4894
Email: pacemedicalinfo@aol.com
www.pacemedicalinc.com

Caution: Federal (USA) law restricts this device to use by or on the order of a physician.

AccuPace™ is a trademark of Pace Medical, Inc.
MICRO-PACE™ is a trademark of Pace Medical, Inc.
Mini Pacer™ is a trademark of Pace Medical, Inc.
Bedside™ is a trademark of Pace Medical, Inc.
Pace Line™ is a trademark of Pace Medical, Inc.
Universal™ is a trademark of Pace Medical, Inc.
Redel® is a registered trademark of Lemo SA.
Pace Line products may use U.S. Patent 6,152,751