

---

## Your Partner...

Pace Medical... over  
30 years of pacing  
experience.

...designs and manu-  
factures Temporary  
Cardiac Pacemakers,  
Pacing Analyzers and  
accessory products  
used in pacing the  
heart, diagnosing the  
performance of im-  
plantable pacing leads  
for permanent pace-  
makers and implant-  
able defibrillators, while  
contributing to better  
patient care.

---

...product performance,  
dedication to quality  
and to the improvement  
of patient care. High  
standards, no compro-  
mise to excellence,  
helping patients one  
day at a time.

---

## ***PACE Medical***

Pace Medical, Inc.  
391 Totten Pond Road  
Waltham, MA 02451  
U.S.A

tel: 781-890-5656  
fax: 781-890-4894  
email: [pacemedicalinfo@aol.com](mailto:pacemedicalinfo@aol.com)  
[www.pacemedicalinc.com](http://www.pacemedicalinc.com)

***PACE Medical***

---

*...Solutions today for  
Your pacing needs...*

*Technology  
Reliability  
Clinical Experience*



---

**Discover the Difference™.**

# Corporate Profile...

Pace Medical is positioned to be the leading medical device company exclusively dedicated to designing and marketing a complete line of high-quality temporary cardiac pacing products to meet the expanding needs of the medical profession worldwide. Our products include both single and dual-chamber temporary cardiac pacemakers, a dual-chamber pacing analyzer, autoclavable and single-use temporary pacemaker patient extension and surgical cables and related accessories.

## Products

The AccuPace™, dual-chamber pacing analyzer, combines the benefits of dual-chamber, multi-mode, multi-parameter pacing with testing features that evaluate, display and store important characteristics of the patient's lead system. The temporary pacing feature allows the physician to perform pre-implant stimulation



AccuPace-  
Pacing Analyzer

studies of the basic pacing parameters and functions on the patient prior to the implantation of a permanent pulse generator. Considerable attention has been paid to the "user-friendly" aspects of operation.

Our single-chamber temporary cardiac pacemakers include the Bedside™ model which attaches to the patient's bedrail and the Mini-Pacer™ which enables the patient to remain ambulatory during treatment.

Both products are used worldwide and have significant clinical use.

**MICRO-PACE-  
Dual-Chamber Temporary  
Pacemaker**



The Company's MICRO-PACE™, dual-chamber temporary cardiac pacemaker, is software programmable and capable of multi-parameter, multi-mode operation. It has the ability to sense and pace at high rates in the DDD mode; thereby, restoring A-V synchrony and improving cardiac output.

The Pace Line™ series of autoclavable and single-use products include— surgical and patient extension cables and adapters— designed to enhance our ability to offer a complete line of highly-competitive temporary pacing products for our customers.

*...Pace Medical, your answer to your pacing needs...*

## **PACE Medical**

Pace Medical, Inc.  
391 Totten Pond Road  
Waltham, MA 02451  
U.S.A

tel: 781-890-5656  
fax: 781-890-4894  
email: [pacemedicalinfo@aol.com](mailto:pacemedicalinfo@aol.com)  
[www.pacemedicalinc.com](http://www.pacemedicalinc.com)